



PRODUCT CATALOGUE
LOW VISION DEVICES

AKRITI OCULOPLASTY LOGISTICS

CORPORATE OFFICE:

8-1-284/OU/232, OU Colony
Shaikpet, Hyderabad - 500 008
Andhra Pradesh, INDIA
Care: +91-995 954 1995,+91-996 318 8447,
Tel: +91-40-32443781, Fax: 91-40-23562729
Email: sales@akriti.co.in
Web : www.akriti.co.in

SALES OFFICE:

6-3-927, Beside Shell Petrol Pump,
Raj Bhavan Road, Somajiguda,
Hyderabad - 500 082
Andhra Pradesh, India
Tele Fax: 91-40-42210407, +91-996 318 8447,
Email: hyderabad@akriti.co.in
Web : www.akriti.co.in

LOW VISION AIDS

HAND HELD MAGNIFIERS




Item Name : Pocket magnifier
Item Code : 0010082
Description : Magnification 6X Field of view 20mm
Working Distance 30mm



Item Name : Hand held magnifier
Item Code : 0010083
Description : Model Magnification Field of View
Working distance L 118 A 3X 20mm 40mm L 118 B
2.5X 30mm 50mm L 118 C 2.5X (5X) 60mm 85mm L
118 D 2.0X 95mm x 65mm 105mm



Item Name : Reading bar
Item Code : 0010084
Description : 2X, Made of acrylic rod



Item Name : Large area illuminated magnifier
Item Code : 0010085
Description : Magnification Field Of view Working
Distance 1.75X (3X) 250mm x250mm 250mm




Item Name : Fresnel wallet magnifier
Item Code : 0010086
Description : Magnification 3X Field of view 65mm x
40mm Working distance 50mm



Item Name : Fresnel book magnifier
Item Code : 0010087
Description : Model L103A L103B Magnification 2.00X
1.75X Field Of view 160mm x 90mm 250mm x 150mm
Working Distance 120mm 250mm



Item Name : Stand magnifier
Item Code : 0010088
Description : Magnification Field Of view Working
Distance 2.5X 30mm 50mm



Item Name : Cutaway stand magnifier
Item Code : 0010097
Description : Magnification Field Of view Working
Distance 6X 22 mm 25mm



Item Name : Wide field magnifier
Item Code : 0010099
Description : Magnification Field Of view Working
Distance 4X 44mm 52mm



Item Name : Pocket magnifier
Item Code : 0010100
Description : Magnification Field Of view Working
Distance 3X 20mm x 20mm 40 mm



Item Name : Hastings triplet loupe
Item Code : 0010102
Description : Model L126A L126B L126C L126D
Magnification 10X 10X 15X 20X Field of View 15mm
20mm 12mm 8mm Working dist. 22mm 22mm 16mm
11mm

PRODUCT DESCRIPTION



Item Name : Folding loupe
Item Code : 0010103
Description : Magnification Field Of view Working Distance 10X 17mm 28mm

SPECTACLE MAGNIFIERS



Item Name : LVA aspheric spectacles
Item Code : 0010104
Description : LVA aspheric spectacle is an useful reading aid for persons with impaired vision. Available powers ($^{\circ}$ D): +16.0, + 20.0, + 24.0, + 28.0



Item Name : LVA prismatic spectacles
Item Code : 0010105
Description : Available powers (D): +6D sph 8D Prism (base inward) +8D sph 10D Prism (base inward) +10D sph 12D Prism (base inward) +12D sph 14D Prism (base inward)




Item Name : Stereo eyeglass magnifier
Item Code : 0010106
Description : Model L 120 A L 120 B Magnification 3X 2.5X Field of View 50mm 50mm Working dist. 340mm 500mm

NONOPTICAL DEVICES



Item Name : Felt tip pen
Item Code : 0010107
Description : Excellent pen for low vision patient to work with great comfort



Item Name : Letter writer guide
Item Code : 0010108
Description : Excellent tool for low vision patients to write in a line work great with dyslexia patients



Item Name : Signature guide
Item Code : 0010109
Description : Excellent tool for low vision patients to make signature on checks write in a line work great with dyslexia patients also

ELECTRONIC DEVICES



Item Name : Hand held illuminated magnifier
Item Code : 0010090
Description : Magnification Field Of view Working Distance 3X 20mm x 20mm 40mm

PRODUCT DESCRIPTION



Item Name : TV filter cum magnifier
Item Code : 0010091
Description : Model L105A L105B L105C
 Magnification 1.5X 1.5X 1.5X Sizes of Lens (Diagonal)
 760mm 690mm 610mm



Item Name : Head loupe
Item Code : 0010092
Description : Model L111A L111B L111C L111D
 Magnification 2.00X 2.75X 3.5X 1.75X Field of View
 120mm 80mm 55mm 140mm Working distance 140mm
 100mm 70mm 180mm



Item Name : Illuminated graticule magnifier
Item Code : 0010093
Description : Magnification, field of view, working
 distance: 8X 2,2 mm, 20mm



Item Name : Measuring microscope
Item Code : 0010094
Description : Model L 119 A L 119 B L 119 C
 Magnification 20X 60X 100X Field of View 9.5mm
 3mm 1.5mm Working dist. 35mm 11mm 5.5mm



Item Name : Pocket microscope
Item Code : 0010095
Description : Magnification Field Of view Working
 Distance 30X 4mm 25mm



Item Name : LVA scope

Item Code : 0010096

Description : A 14' monitor giving crystal clear and sharp images. Superb depth of field, so that all parts of the object remain in sharp focus.