

2011

Ocularist Equipments



Akriti Oculoplasty Logistics

Hyderabad, India

1/1/2011

OCULARISTRY EQUIPMENTS & CONSUMABLES

Code: 0010115: - Impression tray ,

A set of wide range of impression tray made in high quality PMMA, with steady PVC Stem, can be done sterilized with choloro hexidine



Code: 0010049: - Tropicalgin,

Used for taking the imprint of the defects can be used with water/ saline has two color at mixing time pink that turn to yellow working time 1.5 mins



Code: 0010050: - Zelgan

Used for taking the imprint of the defects can be used with water/ saline single color working time 2-3 mins



Code: 0010119: - Carving tool set,

A set of all stainless steel instruments for making prosthesis. can be sterilized by autoclave method



Code: 0010030: - Metal flask large, Price: 150.00, Shipping Cost: 1000.00

Single unit with three pieces, Can be used multiple times



Code: 0010029: - Metal flask Small,

Single unit with two pieces, Can be used multiple times



Code: D-2: - Casting stone,

Come in pack of 1kg and 2 kg Green in Colour



Code: 0010028: - Iris buttons,

Used a clear corneal unit in ocular prosthesis/ scleral shell Flat and Curve units are available



PRODUCT DESCRIPTION

2011

Code: 0010027: - Metal Iris Buttons,

Use in locating the iris plan of the prosthesis, Light in weight non sterile



Code: 0010118: - Metal iris die,

Made up of with high quality stainless steel, single cavity or multiple cavity high shine, high gloss, able to produce the corneal units



Code: 0010028B: - Ocular Black Disc,

A black Ocular Disc used in developing prosthetic eyes



Code: 0010017: - Clear Polymer,

Pack of 300 grams of clear polymer powder , with high molecular weight of Poly methyl metha acrylate



Code: 0010016: - White Polymer Base PMMA,

Pack of 300 grams of white scleral base, with high molecular weight of Poly methyl metha acrylate



Code: 0010055: - Base plate wax, modeling wax,

Flat strips of pink color wax, ideal for the pattern of the prosthesis




Code: 0010053: - Blue inlay wax,

Round stick of 2 incs Hard wax, Cool very quick



Code: 0010057: - Carving wax

Block of white color Hard wax, Cool slowly



PRODUCT DESCRIPTION

2011

Code: 0010056: - Spur wax, wire wax

Round stick of 6.0 incs Blue in color Hard wax, Cool very quick



Code: 0010054: - Sticky wax

Round stick of 2.5 incs Hard wax, Cool very quick pack of 12 pcs



Code: 0010048: - METAL BURR SET

Made of stainless steel Can be sterilized with autoclave Size come from 1.0 mm,1.5mm,2.0mm,2.5mm,3.0mm,3.5mm,4.0mm,4.5mm,5.0mm,5.5mm,6.0mm



Code: 0010117: - Fine Sable Brush

Round ferrule, round or pointed tip. Useful for detail, wash, fills, and thin to thick lines. A pointed round is used for fine detail. A detailer is a pointed round with very short hair.



Code: 0010116: - Dry Color Pigments

wide range of the colors, non toxic, work great with the monopoly syrup wide range



Code: 0010129: - AIRBASE

Hose bench type collector: it has been realized with a special hose,which could become hard till its blockage or even flexible,without the mantainance of position,simply screwing or unscrewing coils handly.Equipped with rubber head and metal base,it could be fixed with bi-adhesiv sticker or with screws,at the bench or on the wall,back to the working bench.Suitable to be used with AIRBOX or MINIBOX EMC



Code: 0010127: - DRILL PLASTER MACHINE PIN 2000

The PIN 2000 Drill Plaster Machine is capable of drilling plaster models or plexiglas and PVC bases for the applications of conical pins in the removable stump technique. The machine is supplied with Laser pointer system that allows a fast and safe alignment of the perforation. This is done by lowering the working surface with both hands. The perforation depth can be regulated by rotating a nut clearly located on the front of the machine.



Code: 0010024: - Emery power,

Ideal for finishing the acrylic products



Code: 0010022: - Emery sheets,

Available no are from 100-1000



Code: 0010060: - Heavy Duty Mandrills,

Made up of stainless steel Heavy Duty Mandrills



Code: 0010131: - HYDRAULIC PRESS

HYDRAULIC PRESS CLOSES UP TO 3 MUFFLES It can reach a max.pressure of 400 bar and produce on the muffles even a max.strength of KG.15.000



Code: 0010058: - Insertion suction cups

Made of high density silicon long lasting Used for inserting the ocular prosthesis



Code: 0010126: - LABO STATION,

LABO STATION is a safety working place cabinet for finishing works with micromotors. It is equipped with several lateral entrances both on the left and on the right side, with a neon light, fast air connector and power supply intake for the dust extractor unit AIR BOX (for simultaneous switch-on/off), wooden air-stock-intake, internal work-top recovered of wood and easily replaceable, a wide sized glass window, a drain trellis for the falling down of heavy powders, and air syringe.



PRODUCT DESCRIPTION

2011

Code: 0010120: - Large pipettes

Pipettes are often called droppers. Use these transfer pipettes to fill lip balm tubes, measure small amounts of fragrance, liquid colorants and more. Once you have these on hand, you will find even more uses. The small dropper has a capacity of 3.5 mL, there are 5 mL in a teaspoon. Markings on the side will show each 1/2 mL increment.




Code: 0010026: - MAXI TORCH EXPERT 27,

Practicle MAXI TORCH EXPERT 27 Piezo electric lighting on,gas lasting: one hour,flame temperature 1300°C rechargeable with normal propane-butane gas



Code: 0010125: - MICROBOX FINISHING CABINET

Cabin for the milling and finishing of precious metals,and for the treatment of hazardous volatile chemical compounds. A pull-out drawer collects residues of the precious alloys being used.The tube supplied with the unit enables to connect an AIRBOX GS or any other aspiraiton system.It can also be used for small milling jobs in dental surge- ries. It is made from plexiglass and all joints are proof bonded.



Code: 0010132: - NATURAL GAS SICURBUNSEN

NATURAL GAS SICURBUNSEN WITH SAFETY THERMOCOUPLE



Code: 0010023: - Pumice Powder,

Fine coarse power for excellent finishing of acrylic work



PRODUCT DESCRIPTION

2011

Code: 0010025: - REFILL CONTAINER,

For Micro Gas Torches filling gas torch



Code: 0010059: - Removal suction cups,

Made of high density silicon rubber long lasting used for removal of prosthesis



Code: 0010020: - Separating sheets

Pack of 100 polyethylene sheets Single use



CODE: 0010130: - T/30 HANDPIECE MM 2,35 COLLET

Laboratory Hand piece,complete with mm.2.35 collet



Code: 0010021: - Vein threads

Pack of 100 meter Single unit



Code: 0010019: - Wax relasing solvent

Pack of 150 ml

Code: 0010128: - WOODEN VACUUM BLOCK

WOOD-WEDGE FOR DUST SUCTION COMPLETE WITH PLEXIGLAS PROTECTION



PRODUCT DESCRIPTION

2011

Code: ERK-1: - Liquisteam-E,

Acryliser / Pressure Polymerization Equipment Model Liquisteam



Code: 0010122: - LUX-1 DENTAL LATHE,

POLISHING LATHE LUX-1,1 SPEED 2.800 RPM 300 WATT WEIGHT 6 KG



Code: 0010123: - LUX-2 DENTAL LATHE

LUX 2 POLISHING LATHE, 2 SPEEDS 1400-2800 RPM,300 WAT



Code: 0010121: - MAP-1 DENTAL LATHE

POLISHING LATHE MAP 1,1 SPEED FOR CONTINUOUS DUTY,500 WATT,2.800 RPM



Code: 0010124: - MINI-MAP-OR DENTAL LATHE

Polishing machine for watch-makers, complete with grinding wheel,spindle and flexible shaft connection. 1 speed (2.800 rpm),200 watt



Code: 0010052: - Liquid Polish

Can be applied over the buffing wheel, to give glass like appearance



Code: 0010051: - Polishing cake

Available in different colors Used for acrylic polishing



PRODUCT DESCRIPTION

2011

Code: 0010042: - Acrylic feathering burr

Small stone acrylic feathering burr Pink in color Can be mounted on any hand piece



Code: 0010040: - Acrylic trimmer

Heavy Duty Blue in color Stone round neck burr with handle can be mounted on any of hand piece.



Code: 0010041: - Acrylic trimmer

Fine trimming burr Flat, Stone round neck burr with handle can be mounted on any of hand piece.



Code: 0010039: - Fast cutting acrylic trimmer

Fast Cutting tool Can be inserted on the head of the lathe machine



Code: 0010045: - Felt Cone Burr

Ideal for working with emery power for acrylic and metal



Code: 0010031: - Silicon fine polisher

Single round silicon polishing burr blue in color, Used for finishing the metal, acrylic products



Code: 0010033: - Silicon fine polisher

Single conical silicon polishing burr blue in color, Used for finishing the metal, acrylic products



Code: 0010032: - Silicon fine polisher

Single conical silicon polishing burr green in color, Used for finishing the metal, acrylic products



Code: 0010034: - Silicon fine polisher

Single round silicon polishing burr green in color, Used for finishing the metal, acrylic products



Code: 0010035: - Silicon fine polisher

Single oval silicon polishing burr blue in color, Used for finishing the metal, acrylic products



PRODUCT DESCRIPTION

2011

Code: 0010036: - Silicon fine polisher,

Single oval silicon polishing burr green in color, Used for finishing the metal, acrylic products



Code: 0010037: - Silicon fine polisher

Single circular burr, 2.5 mm thickness attachable with mandrill Blue in color



Code: 0010038: - Silicon fine polisher

Single circular burr, 2.5 mm thickness attachable with mandrill Green in color



Code: 0010046: - Silicon polisher,

Round silicon polishing burr Gray in color, Made of high density silicon for smoothing of edge of the acrylic shapes



Code: 0010043: - Small Buffing wheels

1 inch diameter, Mounted with small mandrills



Code: 0010047: - Small polishing wheels

Single piece can be mounted with small mandrills can be supplied also with madrils



Code: 0010044: - Soft Buffing wheel

Mounted with mandrill Single piece Can be mounted with any hand piece

