

2011

Punctum Plugs



Akriti Oculoplasty Logistics

Hyderabad, India

1/1/2011

PUNCTAL PLUGS

CODE: 0010112: - DISSOLVABLE COLLAGEN PLUGS (SHORT-TERM)

SHORT-TERM OCCLUSION THERAPY (4-7 DAYS) TEN PACKETS PER BOX, 6 PLUGS PER PACKET. JEWELER'S FORCEPS REQUIRED FOR INSERTION.



CODE: 0010111: - DISSOLVABLE OPAQUE HERRICK LACRIMAL PLUGS ® (MEDIUM)

MEDIUM-TERM OCCLUSION THERAPY (APPROXIMATELY 6 MONTHS) TWO STERILE PLUGS PER BOX, PACKAGED IN A FOIL POUCH TO KEEP PLUGS SOFT, AND PEEL-OFF CHART LABELS FOR ACCURATE DOCUMENTATION, JEWELER'S FORCEPS REQUIRED FOR INSERTION. PURPLE COLOR MAKES THESE PATENTED OPAQUE PLUGS VISIBLE POST-INSERTION WITH TRANSILLUMINATION.



CODE: 0010110: - OPAQUE HERRICK LACRIMAL PLUGS ® NON-DISSOLVABLE (LONG TERM)

TWO INDIVIDUALLY PACKAGED STERILE PLUGS PER BOX, WITH PEEL-OFF CHART LABELS FOR ACCURATE DOCUMENTATION, PRE-MOUNTED ON DISPOSABLE INSERTERS. BLUE COLOR MAKES THESE PATENTED OPAQUE PLUGS UNIQUELY VISIBLE POST-INSERTION WITH TRANSILLUMINATION.

