


2011

Refraction Instrument

AKRITI SUPER SPECIALTY EYE CARE PRODUCTS

Akriti provides all kinds of Pediatric Ophthalmology, Oculoplasty, Cornea & Anterior Segments, Optical, Optometric Products, Low vision, Retina, Glaucoma and Ocularistry Super Specialty Product, Stainless Steel Surgical Sets, Eye Models, Medical Software, Electronic medical record Software (EMR), Hospital Management Systems (HMS), clinic accessories.

Akriti Oculoplasty Logistics
Hyderabad, India



REFRACTION SYSTEM

CODE: DRS3100: - DRS-3100 LITE-CCP HUVITZ DIG REF SYSTEM,

HUVITZ DRS-3100 LITE-CCP THE HUVITZ LITE-CCP INCLUDES THE CDR-3100 DIGITAL REFRACTOR AS WELL AS THE CCP-3100 AUTOMATIC PROJECTOR. THIS PACKAGE PROVIDES THE ABILITY TO ADD DIGITAL REFRACTION TO AN EXISTING EXAM ROOM.



CODE: OH00003: - HUVITZ CCP-3100 AUTO PROJECTOR,

THE HUVITZ CCP-3100 AUTO PROJECTOR IS DESIGNED AS A STAND ALONE AUTOMATED INSTRUMENT THAT ADDS THE CONVENIENCE OF DELIVERING EYE CHARTS TO YOUR PATIENTS VIA A PROGRAMMABLE REMOTE HAND CONTROLLER.



CODE: CPHU02-000: - PROJECTOR HUVITZ CCP-3100 AUTO POC,

THE HUVITZ CCP-3100 AUTO PROJECTOR IS DESIGNED AS A STAND ALONE AUTOMATED INSTRUMENT THAT ADDS THE CONVENIENCE OF DELIVERING EYE CHARTS TO YOUR PATIENTS VIA A PROGRAMMABLE REMOTE HAND CONTROLLER. ISOLATE LETTERS ON A VERTICAL OR HORIZONTAL LINE OR INDIVIDUALLY. THE CCP-3100 IS ALSO DESIGNED TO INTEGRATE WITH HUVITZ'S DRS-3100 DIGITAL REFRACTION SYSTEM.



CODE: BV930: - CHART PROJECTOR -AKRITI

1. PROJECTION DISTANCE:2.0 TO 7.0M
2. PROJECTION ENLARGEMENT RATE:21.5(AT 5M DISTANCE PROJECTION)
3. TILTING ANGLE:5 DEGREE EACH UPWARDS AND DOWNWARDS FROM HORIZONTAL LEVEL
4. LIGHT SOURCE:12V 30W
5. AUTOMATICALLY POWER OFF : 5 MINUTES,10 MINUTES OR 15 MINUTES
6. POWER SUPPLY:10C-220V 50H



CODE: BV935: - CHART PROJECTOR (COOL LIGHT)-AKRITI

AUTO PROJECTOR(COOL LIGHTS) CAN PROVIDE 30 TYPES OF CLEAR ICONS WITH DIFFERENT FUNCTIONS, AND CAN PROGRAM ICON GROUPS PER SPECIFIC NEEDS. IT CAN PRODUCE MORE PRECISE MEASUREMENT RESULTS IN SHORTER TIME BY REMOTE CONTROLLER AND FAST TOGGING ICON, APPLYING RED/GREEN & POLARIZED ICON TEST SHEET, BINOCULAR EYESIGHT AND BINOCULAR BALANCE. USING COOL LIGHTS THIS KIND OF LIGHT CAN BE USED LONG TIME, THE CLEARNESS OF ITS NATURAL LIGHT IS VERY HIGH. MOREOVER, THIS LIGHT IS ENERGY-SAVING. CUSTOMER PROGRAMMABLE ICONS CAN BE DISPLAYED IN PREDEFINED SEQUENCE BY PROGRAMMING. DISPLAY SINGLE ICON DISPLAY SINGLE ICON BY MULTI-FUNCTION MASK PLATE. TILTED PLACEMENT PERFECT IMAGE POSITION CAN BE OBTAINED BY ADJUSTING THE HORIZONTAL AXIS OF PROJECTOR.



CODE: BV932: - CHART PROJECTOR (HALOGEN LAMP) -AKRITI

SPECIFICATION: PROJECTION DISTANCE : 1.5M ~ 6M PROJECTION
MAGNIFICATION : 30 X (AT 5M)



CODE: CDC-4000: - PROJECTOR HUVITZ CDC-4000,

THE HUVITZ CDC-4000 MIRRORED CHART IS THE FUTURE GENERATION OF VISUAL ACUITY TESTING, OVERCOMING THE WEAKNESS OF THE CHART PROJECTOR. IT PROVIDES ACCURATE MEASUREMENTS BY APPLYING CUTTING-EDGE TECHNOLOGY BASED IN AN ERGONOMIC AND COMPACT DESIGN. IT REQUIRES A REFRACTION DISTANCE LESS THAN 4'.



PRODUCT DESCRIPTION

2011

CODE: DRS-3100: - PROJECTOR HUVITZ CDC-4000,

THE HUVITZ CDC-4000 MIRRORED CHART IS THE FUTURE GENERATION OF VISUAL ACUITY TESTING, OVERCOMING THE WEAKNESS OF THE CHART PROJECTOR. IT PROVIDES ACCURATE MEASUREMENTS BY APPLYING CUTTING-EDGE TECHNOLOGY BASED IN AN ERGONOMIC AND COMPACT DESIGN. IT REQUIRES A REFRACTION DISTANCE LESS THAN 4'. THIS ALLOWS A FULL EXAM LANE TO BE SET UP IN A ROOM AS SMALL AS 7' X 7'. THE 30 BUILT-IN CHARTS HAVE BEEN DESIGNED TO MEET YOUR EXACTING NEEDS FOR EXTREME ACCURACY. THE CDC-4000 ALSO DESIGNED TO INTEGRATE WITH MIRAE'S DRS-3100 DIGITAL REFRACTION SYSTEM.



CODE: OH0008: - SIRIUS VISUAL ACUITY SOFTWARE KIT

SOFTWARE AND WIRELESS REMOTE PACKAGE THE SIRIUS VISUAL ACUITY SYSTEM PLACES A WIDE RANGE OF ACUITY TESTS INTO ONE SIMPLE TO USE PRODUCT. SIRIUS RANDOMIZES ALL OPTOTYPES TO PREVENT MEMORIZATION AND PROVIDES A HIGH TECHNOLOGY ADDITION TO YOUR EXAM ROOM. FOR ABOUT THE SAME PRICE AS A MOTORIZED AUTO PROJECTOR THE CHOICE IS CLEAR. CHOOSE THE SIRIUS COMPUTERIZED PROJECTOR FOR SUPERIOR PERFORMANCE. INCLUDE A STANDARD SELECTION COLLECTION OF NATURE VIDEOS, SIMPLE 3D ANIMATIONS FOR CHILDREN, SCENIC VIEWS, PLANET IMAGES, AND OCULAR PATHOLOGY SLIDES.

THE POWERFUL SLIDE SHOW FEATURE ALLOWS YOU TO DEMONSTRATE PROCEDURES, MARKET NEW PRODUCTS, SHOW OCULAR PATHOLOGIES, PRACTICE IMAGES, ETC. THE SIRIUS VISUAL ACUITY SYSTEM CAN BE INSTALLED ON A LAPTOP COMPUTER FOR REMOTE ACUITY TESTING IN SCHOOLS, NURSING FACILITIES, AND HOME VISITS. SIRIUS DUAL MONITOR SOLUTION



CODE: DRS-3100L: - SIRIUS VISUAL ACUITY SYSTEM IMAC

MACINTOSH IMAC TURN-KEY VERSION THE SIRIUS VISUAL ACUITY SYSTEM PLACES A WIDE RANGE OF ACUITY TESTS INTO ONE SIMPLE TO USE PRODUCT. SIRIUS RANDOMIZES ALL OPTOTYPES TO PREVENT MEMORIZATION AND PROVIDES A HIGH TECHNOLOGY ADDITION TO YOUR EXAM ROOM. FOR ABOUT THE SAME PRICE AS A MOTORIZED AUTO PROJECTOR THE CHOICE IS CLEAR. CHOOSE THE SIRIUS COMPUTERIZED PROJECTOR FOR SUPERIOR PERFORMANCE. INCLUDE A STANDARD SELECTION COLLECTION OF NATURE VIDEOS, SIMPLE 3D ANIMATIONS FOR CHILDREN, SCENIC VIEWS, PLANET IMAGES, AND OCULAR PATHOLOGY SLIDES. THE POWERFUL SLIDE SHOW FEATURE ALLOWS YOU TO DEMONSTRATE PROCEDURES, MARKET NEW PRODUCTS, SHOW OCULAR PATHOLOGIES, PRACTICE IMAGES, ETC. THE SIRIUS VISUAL ACUITY SYSTEM CAN BE INSTALLED ON A LAPTOP COMPUTER FOR REMOTE ACUITY TESTING IN SCHOOLS, NURSING FACILITIES, AND HOME VISITS.



CODE: RFCN10-000: - AUTO REFRACTOR, CANON R-F10

THE CANON R-F10 SIMPLIFIES THE STANDARD MEASUREMENT PROCEDURE TO AN UNPRECEDENTED DEGREE - ALL OPERATIONS IN JUST ONE STEP. WITH A PRESS OF THE START BUTTON, THE R-F10 AUTOMATICALLY COMPLETES THE ENTIRE SERIES OF STEPS FOR BOTH EYES: ALIGNMENT, MEASUREMENT, AND PRINTOUT.



CODE: BV610: - AUTO REF-KERATOMETER

SIMPLIFIES THE STANDARD MEASUREMENT PROCEDURE TO AN UNPRECEDENTED DEGREE - ALL OPERATIONS IN JUST ONE STEP. WITH A PRESS OF THE START BUTTON, THE R-F10 AUTOMATICALLY COMPLETES THE ENTIRE SERIES OF STEPS FOR BOTH EYES: ALIGNMENT, MEASUREMENT, AND PRINTOUT.



CODE: RFHU05-000: - AUTO REF/KER WAVEFRONT HUVITZ HRK-7000 ,

THE HRK-7000 WAVEFRONT AUTO REFRACTOR/KERATOMETER IS THE NEWEST GENERATION INSTRUMENT FROM HUVITZ. THIS MODEL COMBINES THE FEATURES OF AN AUTO REFRACTOR, AUTO KERATOMETER AND AN ABERROMETER INTO ONE CONVENIENT-TO-USE INSTRUMENT. THE HRK-7000 IS BASED ON THE HARTMANN-SHACK PRINCIPLE AND FOCUSES ON MEASURING SECOND-ORDER ABERRATIONS, WHICH RESULT IN MYOPIA OR HYPEROPIA AND REGULAR ASTIGMATISM. BECAUSE CORRECTING SECOND-ORDER ABERRATIONS HAS THE HIGHEST IMPACT ON ACUITY, WHICH IS THE EYE'S ABILITY TO DISTINGUISH OBJECT DETAILS AND SHAPE.



CODE: RFHU06-000: - AUTO REF/KER WAVEFRONT HUVITZ HRK-7000A,

THE HRK-7000A WAVEFRONT AUTO REFRACTOR/KERATOMETER IS THE NEWEST GENERATION INSTRUMENT FROM HUVITZ. THIS MODEL COMBINES THE FEATURES OF AN AUTO REFRACTOR, AUTO KERATOMETER AND AN ABERROMETER INTO ONE CONVENIENT-TO-USE INSTRUMENT. THE HRK-7000A IS BASED ON THE HARTMANN-SHACK PRINCIPLE AND FOCUSES ON MEASURING SECOND-ORDER ABERRATIONS, WHICH RESULT IN MYOPIA OR HYPEROPIA AND REGULAR ASTIGMATISM. BECAUSE CORRECTING SECOND-ORDER ABERRATIONS HAS THE HIGHEST IMPACT ON ACUITY, WHICH IS THE EYE'S ABILITY TO DISTINGUISH OBJECT DETAILS AND SHAPE.



CODE: RFR01-000: - AUTO REF/KERATOMETER REICHERT RK600,
 REICHERT RK600 AUTO REFRACTOR / KERATOMETER OVERVIEW THE FULL COLOR LCD SCREEN, ICON-BASED OPERATING SYSTEM, AND FAMILIAR JOYSTICK CONTROLS, ENABLE ANY MEMBER OF YOUR STAFF TO MAKE QUICK AND ACCURATE MEASUREMENTS. MEASUREMENT DATA CAN BE VIEWED ON THE COLOR SCREEN, AND RESULTS CAN BE PRINTED WITH THE FAST AND QUIET, BUILT-IN PRINTER.



CODE: RFNI10-000: - AUTO REF/KERATOMETER RETINOMAX K+ 3,

RETINOMAX K+ 3 BY RIGHT MEDICAL FORMERLY NIKON MOBILE SHARE BETWEEN OFFICES OR BETWEEN ROOMS - THE RETINOMAX K+ 3 IS A FULL FUNCTION HANDHELD AUTOREFRACTOR PLUS A KERATOMETER THAT CAN BE USED ANYWHERE. MEASURE PATIENTS IN ANY POSITION, IN A HOSPITAL BED, IN A WHEELCHAIR OR IN YOUR PRE-TEST ROOM. GREAT FOR SCREENING AT A SCHOOL OR NURSING HOME. BATTERY OPERATION ELIMINATES THE LIMITATION OF BEING PLUGGED IN AND ALLOWS FUNCTIONALITY ANYWHERE.



CODE: RFNI11-000: - AUTO REFRACTOR, RIGHT MED RETINOMAX 3,

RETINOMAX 3 BY RIGHT MEDICAL FORMERLY NIKON MOBILE SHARE BETWEEN OFFICES OR BETWEEN ROOMS - THE RETINOMAX IS A FULL FUNCTION HANDHELD AUTOREFRACTOR THAT CAN BE USED ANYWHERE. MEASURE PATIENTS IN ANY POSITION, IN A HOSPITAL BED, IN A WHEELCHAIR OR IN YOUR PRE-TEST ROOM. GREAT FOR SCREENING AT A SCHOOL OR NURSING HOME. BATTERY OPERATION ELIMINATES THE LIMITATION OF BEING PLUGGED IN AND ALLOWS FUNCTIONALITY ANYWHERE



CODE: HIS-5000: - DIGITAL CAMERA, HIS-5000,
 THE HUVITZ HIS-5000 DIGITAL CAMERA IS DESIGNED TO FIT ON THE HS-5000 OR HS-5500 SLIT LAMPS.



CODE: LEHU03-000: - LENSOMETER, CLM-3100P (A),

THE HUVITZ CLM-3100P AUTO-LENSMETER EASILY MEASURES PROGRESSIVE MULTI-FOCAL LENSES. THE DATA IS REPRESENTED IN A SIMPLE TO UNDERSTAND GRAPHIC DISPLAY.



PRODUCT DESCRIPTION

2011

CODE: BV560: - AUTO LENSMETER,

SPECIFICATIONS SPHERE: MEASURING RANGE-25~+25D UNIT0.01/0.12/0.25D
CYLINDER: MEASURING RANGE-9.99~+9.99D UNIT0.01/0.12/0.25D AXIS:
MEASURING RANGE 0°~180°UNIT 1° ADD: MEASURING RANGE0~+9.99D
UNIT0.01/0.12/0.25D PRISM: MEASURING RANGE0~0.99Δ UNIT0.01/0.12/0.25Δ
PD: MEASURING RANGE 40~86MM UNIT0.5MM DIAMETER: Ø5~120MM DISPLAY:
LCD MONITOR (5.7INCHES) POWER SUPPLY: AC100W 50/60HZ 40VA



CODE: 93.010: - OPERATING MICROSCOPES,

AB SERIES OF OPERATING MICROSCOPES OF QUALITY RELIABLE PRODUCT THAT
WILL GIVE YOU GOOD MAGNIFICATIONS, ILLUMINATION AND DEPTH. THE
AKRITIM MICROSCOPE PROVIDES SHARP IMAGE RESOLUTION WITH HIGH
INTENSITY COAXIAL ILLUMINATION FOR A WIDE RANGE OF SURGICAL
PROCEDURES INCLUDING EYE, ENT, NEURO-SURGERY, GYNECOLOGY, PLASTIC
SURGERY AS WELL AS EXTRA CAPSULAR CATARACT EXTRACTION AND LENS
IMPLANTATION BESIDES ADVANCED DIAGNOSTIC APPLICATIONS.



CODE: 31200: - OPERATIONG MICROSCOPE (31200),

FIVE-STEP MAGNIFICATIONS FOR MAIN MICROSCOPE, THREE-STEP
MAGNIFICATIONS FOR ASSISTANT MICROSCOPE.



CODE: SP-100: - PACHYETER TOMEY,

A SIMPLE HANDY PACHYETER THE SP-100 IS WHAT YOU WANT AND TOMEY
DELIVERS. THE TOMEY SP-100 HANDY PACHYETER IS ONLY 530G (19 OZ.) AND IS
THE BEST EQUIPMENT TO MEASURE CENTRAL CORNEAL THICKNESS.



CODE: FEOP02-000: - PERIMETER, AP2000 AUTOMATED,

AUTOMATED PERIMETER - AP 2000 EYE TRACKER SYSTEM FOR RAPID EXAMS &
ACCURATE RESULTS.



PRODUCT DESCRIPTION

2011

CODE: TPVE08-000: - RESEEVIT PULSAR TOMOGRAPHY PACKAGE,

THE RESEEVIT PULSAR TOMOGRAPHY PACKAGE COMBINES THE TECHNOLOGIES OF A PLACIDO DISK AND A 3D ROTATING SCHEIMPFLUG CAMERA TO PROVIDE UNPARALLELED DETAIL OF THE CORNEA AND ANTERIOR CHAMBER. DATA FROM 25 SCHEIMPFLUG FRAMES AS WELL AS A 22 RING PLACIDO DISK ARE ACQUIRED IN LESS THAN ONE SECOND AND COMPUTER ANALYZED. BY COMBINING THE INFORMATION FROM THE TWO METHODOLOGIES EXTREME ACCURACY IS ACHIEVED.



CODE: 23.060: - SLIT LAMP ECONOMY,

SLIT LAMPS UTILIZE HIGH QUALITY OPTIC, PROVIDING CLARITY AND DETAIL THAT WILL SURPRISE EVEN THE MOST DISCRIMINATING OPERATOR. THE XYZ JOYSTICKS, CONVENIENT CONTROLS, AND ERGONOMIC DESIGN, COMBINED WITH TOP QUALITY MATERIALS AND METICULOUS ASSEMBLY, ENSURE YEAR OF COMFORTABLE, ENJOYABLE USE. THE OUTSTANDING OPTICAL PERFORMANCE, FEATURES AND QUALITY MAKE THE AKRITI SURGICAL SERIES COMPARABLE TO THE WORLD'S ELITE LAMPS.



CODE: HS-5000: - SLIT LAMP, HUVITZ HS-5000 ,

THE HUVITZ HS-5000 DELIVERS OUTSTANDING OPTICAL QUALITY AT AN AFFORDABLE PRICE POINT. THE NEWEST GENERATION INSTRUMENT FROM HUVITZ IS DESIGNED TO FILL A MARKET GAP.



CODE: BV235: - CROSS CYLINDER,

SPECIFICATION THE CROSS-CYLINDER LENS FOR REFINING ASTIGMATISM CORRECTIONS. --FULL DIAMETER GLASS LENS --LONG HANDLE --SILVER METAL RIM --WOODEN CASE --DIOPTERS: -0.25, -0.50, -0.75,-1.00



CODE: 100070: - DISTOMETER SET,

MEASURES THE VERTEX DISTANCE BETWEEN CORNEA AND TRIAL LENS OR EYEGLASS LEVEL.



PRODUCT DESCRIPTION

2011

CODE: BV561: - MANUAL LENSMETER (EXTERIOR READING),



CODE: BV562: - MANUAL LENSMETER (INTERNAL READING)



CODE: 32.010: - KERATOMETER

KERATOMETER IS THE WORLD WIDE STANDARD FOR MANUAL MEASUREMENT OF CORNEAL CURVATURES.



CODE: BV530: - KERATOMETER

KERATOMETER WAS DEVELOPED UNDER THE GUIDANCE OF OPTOMETRY CENTER. KERATOMETER WAS DEVELOPED UNDER THE GUIDANCE OF OPTOMETRY CENTER. THIS INSTRUMENT CAN GIVE A PRECISE EXAMINATION OF THE RADIUS OF CURVATURE OF THE ANTERIOR CORNEAL SURFACE. THIS MEASUREMENT IS UTILIZED TO FIT CONTACT LENSES AND TO MONITOR CORNEAL CHANGES PRODUCED THROUGH THE WEAR OF CONTACT LENSES. KERATOMETER IS OCCASIONALLY USED TO ASSIST IN THE RECOGNITION OF CERTAIN COMEAL ABNORMALITIES AND TO CHECK THE RADII OF CURAVATURE OF BOTH HARD AND SOFT CONTACT LENSES. THE PROGRAM CAN BE OPERATED WITH CONSIDERABLE EASE.



CODE: BV156: - PD METER ,

SPECIFICATION: 1.DATE RANGE : BINOCULAR PD:46MM-82MM MONOCULAR PD: 23MM-41MM 2.DISTANCE RANGE: 300-INFINITY~ 3.INDICATION ERROR: =0.5MM OR <0.5MM 4.ASYMMETRY ERROR:=0.5MM OR <0.5MM 5.ROUNDING ERROR: =0.5MM OR <0.5MM 6.POWER SOURCE: DC 6V (4PIECES AA BATTERIES)



CODE: BV612: - PHOROPTER,



PRODUCT DESCRIPTION

2011

CODE: BV123: - TRIAL FRAME,

MAIN TECHNICAL INDEXES *ADJUSTMENT OF PUPILLARY DISTANCE: TWO EYES' PUPILLARY DISTANCE IS 54 MM~70 MM, AND ADJUSTMENT OF EACH GRADE REQUIRES CHANGE OF 2 MM. *AXIS GRADUATION OF CALIBRATED DIAL: LEFT AND RIGHT CALIBRATED DIALS ARE IDENTICAL, WHERE AXIS GRADUATION LINE RUNS CLOCKWISE ALONG AXIS GRADUATION LINE OF LENS FRAME, 165°, ~0°, ~180° AND ~25°, WITH DISTANCE OF GRADUATION LINE BEING 5°. *INNER PORE DIAMETER OF LENS FRAME: Φ 34 MM *NUMBER OF LENS THAT CAN BE INSERTED INTO LEFT AND RIGHT LENS FRAMES AT THE SAME TIME: 4 PIECES *DEGREE THAT LENS ROTATES AROUND OPTIC AXLE INSIDE LENS FRAME: 360°



CODE: BV128 : - TRIAL FRAME,

SPECIFICATIONS & OPERATION INSTRUCTIONS: USAGES: IT IS SUITABLE FOR SPECTACLES SHOPS, DEPARTMENT OF OPHTHALMOLOGY IN HOSPITALS, SCHOOLS AND ENTERPRISES TO MEASURE EYESIGHT.



CODE: BV202: - TRIAL LENS SET,

SPECIFICATION: CHARACTER: LENSES WITH PLASTIC AND METAL COMBI, RIM. NOT INCLUDING TRIAL FRAME. SPHERES : 40 PAIRS EACH OF CONCAVE AND CONVEX 0.25D TO 6.00D IN 0.25 STEPS 6.50D TO 10.00D IN 0.50 STEPS 11.00D TO 15.00D IN 1.00 STEPS 16.00D TO 20.00D IN 2.00 STEPS CYLINDERS : 20PAIRS EACH OF CONCAVE AND CONVEX 0.25D TO 6.00D IN 0.25 STEPS 4.50D TO 6.00D IN 0.50 STEPS PRISMS : 12 PIECES 0.5(2) 1.0 TO 10.0 IN 1.0 STEPS ACCESSORIES : 14 PIECES RF GF BL PL(2) MR(2) FL CL(2)PH (2) SS PF CROSSCYLINDER X 0.25 X0.50



CODE: BV205: - TRIAL LENS SET,

COMPOSITION: CHARACTER: LENSES WITH COLOR PLASTIC RIM, LEATHER CASE AND WOODEN TRAY. NOT INCLUDING TRIAL FRAME. SPHERES : 40 PAIRS EACH OF CONCAVE AND CONVEX 0.25D TO 6.00D IN 0.25 STEPS 6.50D TO 10.00D IN 0.50 STEPS 11.00D TO 15.00D IN 1.00 STEPS 16.00D TO 20.00D IN 2.00 STEPS CYLINDERS : 20PAIRS EACH OF CONCAVE AND CONVEX 0.25D TO 6.00D IN 0.25 STEPS 4.50D TO 6.00D IN 0.50 STEPS PRISMS : 12 PIECES 0.5(2) 1.0 TO 10.0 IN 1.0 STEPS ACCESSORIES : 14 PIECES RF GF BL PL(2) MR(2) FL CL(2)PH (2) SS PF CROSSCYLINDER X 0.25 X0.50



PRODUCT DESCRIPTION

DIAGNOSTIC DEVICES

CODE: 11772: - COAXIAL OPHTHALMOSCOPE,

THE BEST OF THE BEST IN TRADITIONAL OPHTHALMOSCOPES COAXIAL VISION SYSTEM, SEALED OPTICS, AND HALOGEN HPX™ ILLUMINATION FOR EXCELLENT VIEWS OF THE FUNDUS.



CODE: 12500: - OPHTHALMOSCOPES / RETINOSCOPES BINOCULAR INDIRECT ,

THE TOP OF A SUPERIOR LINE OF OPHTHALMIC INSTRUMENTS THE WELCH ALLYN BIO IS DESIGNED TO ACHIEVE PRECISE BINOCULAR VIEWING AND STEREOPSIS EASIER THAN EVER BEFORE. ITS OBSERVATION PATHS CAN BE ADJUSTED FROM 48 MM TO 74 MM TO HANDLE THE WIDEST AND NARROWEST INTERPUPILLARY DISTANCES.



CODE: 12860: - OPHTHALMOSCOPES / RETINOSCOPES POCKETSCOPE™ OPHTHA,

WELCH ALLYN QUALITY IN A CONVENIENT POCKET SIZE THE POCKETSCOPE™ FITS INTO YOUR POCKET! EASY TO USE, LIGHTWEIGHT, AND PERFECT FOR GENERAL AND SPECIALIST USE. TAKE YOUR POCKETSCOPE™ WITH YOU EVERYWHERE YOU GO THROUGHOUT THE DAY.



CODE: 11810: - PANOPTIC™ OPHTHALMOSCOPE

INNOVATIVE OPTICS FOR PANORAMIC VIEWS OF THE FUNDUS WELCH ALLYN'S PANOPTIC™ OPHTHALMOSCOPE IS A REVOLUTIONARY NEW CONCEPT FOR VISUALIZING THE FUNDUS. DEVELOPED BY WELCH ALLYN®, THE ORIGINATOR OF THE WORLD'S FIRST HAND-HELD DIRECT-ILLUMINATING OPHTHALMOSCOPE, THE PANOPTIC OPHTHALMOSCOPE INCORPORATES AN INNOVATIVE AND PATENTED NEW OPTICAL DESIGN CALLED AXIAL POINTSOURCE™ OPTICS.



CODE: 12850: - POCKET JUNIOR OPHTHALMOSCOPE

LIGHTWEIGHT AND COMPACT. HALOGEN ILLUMINATION PROVIDES BRIGHTER, WHITER, LONGER LASTING LIGHT. SEALED OPTICE KEEP INSTRUMENT CLEAN AND DUST FREE . APERTURE DIAL OFFERS SIX SELECTIONS. GLASS CONDENSING LENS. INCLUDES RED-FREE AND COBALT BLUE FILTER.



CODE: 18342: - RETINOSCOPES,

SUPERIOR ILLUMINATION, THE RIGHT DESIGN FOR ACCURATE REFRACTION. THE 3.5V RETINOSCOPE WHICH PROVIDES SUPERIOR ILLUMINATION, IS THE RIGHT DESIGN FOR ACCURATE REFRACTION. IT HAS THE UNIQUE QUALITY OF ONE-HAND PERATION AND THE COMFORT OF A RUBBER BROW REST.

