


## AIR INJECTION &amp; IRRIGATING CANNULA




**Item Name :** IRRIGATING CANNULA  
STRAIGHT BLUNT TIP

**Item Code :** MM049

[Enquire.....](#) [Read More.....](#)

**Description :**

20G/ 22/23/25/27/ 30G x 25 mm Length  
MM049-S Set of 6



**Item Name :** ANTERIOR CHAMBER  
CANNULA (RYCROFT)

**Item Code :** MM051

[Enquire.....](#) [Read More.....](#)

**Description :**

Angled 45° overall length 22 mm

MM051-21: 21 G


MM051-23 : 23 G

MM051-25 : 25 G

MM051-26 : 26 G

MM051-27 : 27 G

MM051-30 : 30 G



**Item Name :** BISHOP HARMON  
CANNULA

**Item Code :** MM055

[Enquire.....](#) [Read More.....](#)


**Description :**

40° Angled, spatulated, 8 mm from end,

Length 25 mm

MM055-16 16

MM055-19 19 G



**Item Name :** FLATTENED ANTERIOR  
CHAMBER CANNULA (BRACKEN)

**Item Code :** MM057

[Enquire.....](#) [Read More.....](#)

**Description :**

Curved 90° flattened bevel tip 19 G  
Length 23 mm


**Item Name :** ANTERIOR CHAMBER  
CANNULA (MCINTRYRE)

**Item Code :** MM061

[Enquire.....](#) [Read More.....](#)

**Description :**

MM061Set : Set of 5 Cannulae  
MM061 26(A) : 26 G Angled  
MM061 30(A) : 30 G Angled  
MM061-25(A) : 25 G Angled  
MM061-27(A) : 27 G Angled  
MM061-23(A) : 23 G Angled




**Item Name :** ANTERIOR CHAMBER  
CANNULA (KNOLLS)

**Item Code :** MM063

[Enquire.....](#) [Read More.....](#)

**Description :**

2 mm Angled Tip, Length 19 mm  
MM063 30, 30 G  
MM063 23, 23 G  
MM063 25, 25 G  
MM063-Set Set of 3 Cannulae




**Item Name :** ANTERIOR CHAMBER  
MAINTAINER (LEWICKY)

**Item Code :** MM065

[Enquire.....](#) [Read More.....](#)

**Description :**


Silicone tubing attached to 20G x 3.50 mm  
cannula  
unique annular rings for self retaining



**Item Name :** VISCO FLOW CANNULA

**Item Code :** MM1056

[Enquire.....](#) [Read More.....](#)




**Description :**

For easier injection of viscoelastic to all meridians Dome shape provides safe insertion  
27G x 24 mm formed

**Item Name :** ANTERIOR CHAMBER MAINTAINER (BLUMENTHAL)

**Item Code :** MM1065

[Enquire.....](#) [Read More.....](#)



**Description :**

20 cms Silicone tubing attached to reinforced  
23G x 5.00 mm cannula with Atkinson point  
Slightly spatulated for self retaining

**Item Name :** A/C MAINTAINER BULBOUS TIP

**Item Code :** MM1067

[Enquire.....](#) [Read More.....](#)



**Description :**

Bulbous tip allows easy insertion and  
minimum tissue disturbance  
20 G x 4 mm from tip

## ANESTHESIA NEEDLES

**Item Name :** RETROBULBAR NEEDLE (ATKINSON)

**Item Code :** MM001

[Enquire.....](#) [Read More.....](#)



**Description :**

MM 001-23 : 23G x 38 mm  
MM 001-25 : 25G x 38 mm  
MM 001-27 : 27G x 38 mm


**Item Name :** RETROBULBAR NEEDLE  
(SHARP POINT)

**Item Code :** MM003

[Enquire.....](#) [Read More.....](#)

**Description :**

MM003-23 : 23G x 38 mm  
MMJ003-25 : 25G x 38 mm  
MM003B-25 : 25G x 50 mm




**Item Name :** PERIBULBAR NEEDLE  
(ATKINSON)

**Item Code :** MM005

[Enquire.....](#) [Read More.....](#)

**Description :**

MM005-23 : 23G x 22 mm  
MM005-25 : 25G x 22 mm  
MM005-27 : 27G x 22 mm




**Item Name :** SUBTENON'S ANESTHESIA

**Item Code :** MM009

[Enquire.....](#) [Read More.....](#)

**Description :**

Curved flattened tip 19G x 25 mm




**Item Name :** SUBTENON'S ANESTHESIA  
(TRI PORT)


**Item Code :** MM011

[Enquire.....](#) [Read More.....](#)

**Description :**

Curved : MM011C-19 : 19G x 25 mm  
Curved : MM011C-21 : 21G x 25 mm  
Angled : MM011A-19 19G x 25 mm  
Angled : MM011A-21 21G x 25 mm






**Item Name :** SUBTENON'S ANESTHESIA (PARA) NEEDLE

**Item Code :** MM013

[Enquire.....](#) [Read More.....](#)

**Description :** 21G, 0.30 mm side port



**Item Name :** RETROBULBAR NEEDLE (CURVED) (ATKINSON)

**Item Code :** MM1003

[Enquire.....](#) [Read More.....](#)


**Description :**

Design facilitates to course around the globe and parallel to the optic nerve. 11 mm bend.

MM1003-25 : 25G x 35 mm

MM1003-23 : 23G x 35 mm

### BIMANUAL ASPIRATION



**Item Name :** ASPIRATING HAND PIECE


**Item Code :** MM351

[Enquire.....](#) [Read More.....](#)

**Description :**

0.40 mm port with polisher attached Cannula

MM351-T : Titanium



**Item Name :** IRRIGATING HAND PIECE


**Item Code :** MM353

[Enquire.....](#) [Read More.....](#)

**Description :**

Attached Cannula has cross 0.50 mm irrigating ports

MM353-T : Titanium




**Item Name :** ASPIRATING CANNULA WITH POLISHER

**Item Code :** MM357

**Description :**

0.40 mm port with polisher



**Item Name :** IRRIGATING CANNULA


**Item Code :** MM359

[Enquire.....](#) [Read More.....](#)

**Description :**

Cannula has crossed 0.50 mm irrigating ports

## CAPSULE POLISHERS




**Item Name :** OLIVE TIP POLISHER  
(PEARCE) (TOP PORT)

**Item Code :** MM1157

[Enquire.....](#) [Read More.....](#)

**Description :**

Olive shaped tips with 0.30 mm top port micro  
sand blasted  
25G x 25 mm




**Item Name :** CAPSULE POLISHER (ANTS)

**Item Code :** MM155

[Enquire.....](#) [Read More.....](#)

**Description :**

Sand blasted 0.3 mm top irrigating port, 22  
mm length  
MM155-S : 23G Straight  
MM155-A : 23G Angled




**Item Name :** CAPSULE POLISHER  
(JENSON)

**Item Code :** MM157

[Enquire.....](#) [Read More.....](#)

**Description :**

All around sand blasted bulbous tip, Length 28 mm, gently curved  
MM157-27 : 27 G  
MM157-25 : 25 G  
MM158-27 : Bottom sand blasted




**Item Name :** POSTERIOR CAPSULE POLISHER (FREEMAN)

**Item Code :** MM159

[Enquire.....](#) [Read More.....](#)

**Description :**

Irrigating port on both sides  
19G, Length 32 mm  
Angled, Sand blasted




**Item Name :** POSTERIOR CAPSULE POLISHER (SIMCOE)

**Item Code :** MM161

[Enquire.....](#) [Read More.....](#)

**Description :**

Reverse 0.30 mm aperture, 16 mm long  
23G Curved, Sand blasted




**Item Name :** POSTERIOR CAPSULE POLISHER(SIMCOE)

**Item Code :** MM163

[Enquire.....](#) [Read More.....](#)

**Description :**

Top port 0.30 mm aperture 16 mm long  
23G curved, Sand blasted




**Item Name :** ANTERIOR CAPSULE SCRUBBER (SIMCOE)

**Item Code :** MM165

[Enquire.....](#) [Read More.....](#)

**Description :**

Tip is Angled 45°, and roughened slightly  
23G, Length 19 mm




**Item Name :** POSTERIOR CAPSULE SCRUBBER (SIMCOE)

**Item Code :** MM167

[Enquire.....](#) [Read More.....](#)

**Description :**

Tip is Angled 45°, and roughened slightly  
23G, Length 19 mm




**Item Name :** CAPSULE POLISHER (KRATZ)

**Item Code :** MM169

[Enquire.....](#) [Read More.....](#)

**Description :**

0.30 mm top port, Angled, Blasted 3 mm from end  
15-169-21 : 21G x 22 mm  
15-169-23 : 23G x 22 mm



**Item Name :** COLLAR BUTTON (GREATHER)

**Item Code :** MM171

[Enquire.....](#) [Read More.....](#)

**Description :**

Micro Iris retractor sand blasted posterior surface for capsule polishing side irrigation,  
Length 25  
MM171-S : 23G Straight  
MM171-A : 23G Angled


**Item Name :** FLEXIBLE TIP POLISHER  
(ENGELS)

**Item Code :** MM173

[Enquire.....](#) [Read More.....](#)

**Description :**

Curved, 23 mm sleeve, length 22 mm  
MM173 : 27G




**Item Name :** FLEXIBLE TIP POLISHER

**Item Code :** MM175

[Enquire.....](#) [Read More.....](#)

**Description :**

Anterior Capsule Polisher  
Angled : 45°  
Bend : 9 mm  
Sleeve Length : 12 mm  
MM175 : 23G x 22 mm



## CORTEX ASPIRATION


**Item Name :** CORTEX ASPIRATING  
CANNULAE

**Item Code :** MM079

[Enquire.....](#) [Read More.....](#)

**Description :**

(MCINTYRE BINKHORST)  
Blunt Tip, Front opening, Hooked, 22 mm  
length  
MM079-R : 26G Right  
MM079-L : 26G Left



**Item Name :** CORTEX ASPIRATOR  
(ALPERN)

**Item Code :** MM081

[Enquire.....](#) [Read More.....](#)


**Description :**

Angled 45°, 5 mm from end

Hook 1.50 mm x 1.75 mm

MM081-R : 25G x 22 mm (Right)

MM081-L : 25G x 22 mm (Left)



**Item Name :** ASPIRATING CANNULA  
(GILLS)

**Item Code :** MM085

[Enquire.....](#) [Read More.....](#)


**Description :**

MM085-5 : Angled to tip 5 mm

MM085-6 : Angled to tip 6 mm

MM085-8 : Angled to tip 8 mm

MM085-10 : Angled to tip 10 mm



**Item Name :** ASPIRATION CANNULA  
(HEALON)

**Item Code :** MM087

[Enquire.....](#) [Read More.....](#)


**Description :**

Angled 10 mm from tip, Length 32 mm

MM087-19 : 19 G

MM087-21 : 21 G

MM087-22 : 22 G



**Item Name :** CORTEX EXTRACTOR  
CANNULA (SIMCOE)

**Item Code :** MM091


[Enquire.....](#) [Read More.....](#)

**Description :**


Smooth closed tip with top aspirating port into curve

MM091-13 : 13 mm long (23 G)

MM091-21 : 21 mm long (23 G)



MM091-13 : 13 mm long (21 G)  
MM091-21 : 21 mm long (21 G)



**Item Name :** CORTEX EXTRACTOR  
CANNULA (SIMCOE-WELSH)

**Item Code :** MM093


[Enquire.....](#) [Read More.....](#)

**Description :**

Angled 10 mm from tip, 0.3 mm, aspirating  
port

MM093-23 : 23 G

MM093-21 : 21 G



**Item Name :** CORTEX EXTRACTOR  
CANNULA (WELSH)

**Item Code :** MM095

[Enquire.....](#) [Read More.....](#)

**Description :**


Flat Olive Tip length 25 mm

MM095-S : 25G Straight

MM095-L : 25G Left

MM095-R : 25G Right

MM095 Set Set of 3 Cannulae



**Item Name :** ALPHA CHROMOSTRYPSIN  
ASPIRATING CANNULA

**Item Code :** MM097

[Enquire.....](#) [Read More.....](#)

**Description :**

(TROUTMAN)

Smooth bulbous tip length 26 mm, 25 G

**Item Name :** U SHAPED ASPIRATING CANNULA

**Item Code :** MM099

[Enquire.....](#) [Read More.....](#)



**Description :**

Hooked 1.50 mm x 1.75 mm wide  
Length : 22 mm  
MM099-21 : 21 G  
MM099-23 : 23 G  
MM099-25 : 25G

**Item Name :** CORTEX EXTRACTOR (CHARLEUX)

**Item Code :** MM1098

[Enquire.....](#) [Read More.....](#)



**Description :**


Bulbous tip allows easy insertion and minimum tissue disturbance.

## CYCLODIALYSIS

**Item Name :** CYCLODIALYSIS CANNULA (CASTROVIEJO)

**Item Code :** MM215

[Enquire.....](#) [Read More.....](#)




**Description :**

15 mm Angled tip, Opening at top  
Semi-sharp sides, 19 Gauge

**Item Name :** CYCLODIALYSIS CANNULA (KELMAN)

**Item Code :** MM219

[Enquire.....](#) [Read More.....](#)



**Description :**

Spatulated tip left side opening  
19G MM219-R Right-opening  
MM219-L Left-opening


**Item Name :** CYCLODIALYSIS CANNULA (SHEETS)

**Item Code :** MM221

[Enquire.....](#) [Read More.....](#)

**Description :**

Flattened 12 mm tip  
Irrigating port on both sides and in front  
Overall length 32 mm, 19 Gauge




**Item Name :** CYCLODIALYSIS CANNULA (GANS)

**Item Code :** MM223

[Enquire.....](#) [Read More.....](#)

**Description :**

Opening on left side, 12 mm Angled tip, 19 Gauge




**Item Name :** CYCLODIALYSIS CANNULA (GANS)

**Item Code :** MM225

[Enquire.....](#) [Read More.....](#)

**Description :**

Both sides opening, 12 mm Angled tip, 19 Gauge




**Item Name :** CYCLODIALYSIS CANNULA (RANDOLPH)

**Item Code :** MM227

[Enquire.....](#) [Read More.....](#)

**Description :**

Opening in front for air injection 19 Gauge  
12 mm Angled tip  
Overall length 29 mm



## CYSTOTOMES

**Item Name :** CYSTOTOME (STRAIGHT)

**Item Code :** MM015D

[Enquire.....](#) [Read More.....](#)

**Description :**




Sterile Disposable 16mm long  
MM015-30 : 30G  
MM015-27 : 27G  
MM015-25 : 25G

**Item Name :** IRRIGATING CYSTOTOME (FORMED)

**Item Code :** MM017D

[Enquire.....](#) [Read More.....](#)

**Description :**




Sterile Disposable  
MM017-30 : 30G x 16 mm  
MM017-26 : 26G x 16 mm  
MM017-27 : 27G x 16 mm  
MM017-25 : 25G x 16 mm

**Item Name :** CYSTOTOME (ANGLED)

**Item Code :** MM019D

[Enquire.....](#) [Read More.....](#)

**Description :**



Sterile Disposable 45°  
Tip rotated 45°  
MM019-25 : 25G  
MM019-27 : 27G  
MM019-30 : 30G x 16 mm

**Item Name :** IRRIGATING CYSTOTOME (REVERSE)


**Item Code :** MM021

[Enquire.....](#) [Read More.....](#)

**Description :**



Sterile Disposable  
Tip bent towards bevel  
25G x 22 mm




**Item Name :** SHARP IRRIGATING CYSTOTOME (KELMAN)

**Item Code :** MM023D

[Enquire.....](#) [Read More.....](#)

**Description :**

Curved sharp cutting edge 25G x 32 mm



**Item Name :** IRRIGATING CYSTOTOME (BLUMENTHAL)


**Item Code :** MM025D

[Enquire.....](#) [Read More.....](#)

**Description :**

Sterile Disposable  
Formed 8.5 mm from bend to tip 25G x 32 mm

## HYDRODISSECTION & HYDRODELINEATION




**Item Name :** NUCLEOUS HYDRODISSECTOR (RAININ)

**Item Code :** MM069

[Enquire.....](#) [Read More.....](#)

**Description :**

Angled 45°, Bend 7 mm, Length 22 mm  
Spatulated on horizontal / vertical plane  
MM069-30 : 30 G  
MM069-27 : 27 G



**Item Name :** NUCLEOUS HYDRODISSECTOR (MCINTYRE) SPATULATED

**Item Code :** MM071

[Enquire.....](#) [Read More.....](#)

**Description :**

Angled 45°, Bend 11 mm, Length 22 mm  
Flattened on horiz. plane

MM071-23 : 23 G  
 MM071-27 : 27 G  
 MM071-30 : 30 G


**Item Name :** J SHAPED MICRO  
 HYDRODISSECTOR (PEARCE)

**Item Code :** MM073

[Enquire.....](#) [Read More.....](#)

**Description :**

Hook 1.75 mm wide 2.25 mm  
 Long length : 22 mm  
 MM073-30 : 30 G  
 MM073-27 : 27 G




**Item Name :** NUCLEOUS HYDROLYSIS  
 NEEDLE (P.KOCH)

**Item Code :** MM075

[Enquire.....](#) [Read More.....](#)

**Description :**

Short bevel tip penetrates the nucleus for  
 manual sculpting to improve lamellar  
 separation  
 Flattened & Smooth insertion under anterior  
 capsule  
 Curved : 25 G  
 Length : 22 mm




**Item Name :** HYDRODISSECTION  
 CANNULA (CHANG)

**Item Code :** MM080

[Enquire.....](#) [Read More.....](#)

**Description :**


27G, 90° angle, flat tip with a beveled  
 opening  
 Flat tip for easy insertion under the capsular  
 rim, provide broad stream of fluid  
 Beveled tip used to rotate the nucleus to ensure  
 that it is completely separated



**Item Name :** MICRO J SHAPE  
HYDRODISSECTOR (PEARCE)

**Item Code :** MM1073

[Enquire.....](#) [Read More.....](#)



**Description :**

Micro size spatulated tip for easy entry  
through wound under capsule Hook : 2.00 mm  
x 2.50 mm

MM1073-27 27G x 22 mm

MM1073-25 : 25G x 22 mm

**Item Name :** NUCLEOUS HYDRODISSECTOR  
(PEARCE)

**Item Code :** MM1075

[Enquire.....](#) [Read More.....](#)




**Description :**

Close tip to prevent clogging while entry.  
Three Irrigating ports for improved fluid  
distribution and lamellar separation.  
25G x 22 mm

**Item Name :** NUCLEOUS  
HYDRODELINEATOR (FEASTER)

**Item Code :** MM1077

[Enquire.....](#) [Read More.....](#)



**Description :**

Closed wedge shaped tip allows easy penetration  
between nucleous & epinucleous.


Two irrigating ports for fluid distribution in both  
the planes

25 G x 22 mm

**Item Name :** TAPERED HYDRODELINEATOR  
(BLUMENTHAL)


**Item Code :** MM1078

[Enquire.....](#) [Read More.....](#)



**Description :**

Tapered tip allows more flow of fluid through a smaller tip.  
Easy insertion under capsule for hydrodissection or into nucleus for hydrodelineation.  
25G x 22 mm



**Item Name :** NUCLEUS HYDRODISSECTOR

**Item Code :** MM1079

[Enquire.....](#) [Read More.....](#)

**Description :**

Horizontally Flattened Tip Size  
27 G x 22 mm

## IOL & NUCLEUS MANIPULATION



**Item Name :** IOL MANIPULATOR  
(JAFFE/MALTZMAN)

**Item Code :** 127-01

[Enquire.....](#) [Read More.....](#)

**Description :**

Fits 0.2 mm diameter holes or larger, Angled 10 mm



**Item Name :** IOL MANIPULATOR (LESTER)

**Item Code :** 127-02

[Enquire.....](#) [Read More.....](#)

**Description :**

Hour Glass' shaped blunt tip, angled 10 mm



**Item Name :** IOL MANIPULATOR (SINSKEY)

**Item Code :** 127-03

[Enquire.....](#) [Read More.....](#)

**Description :**

Fits 0.2 mm diameter holes or larger, Straight

**Item Name :** IOL MANIPULATOR (SINSKEY)

**Item Code :** 127-04

[Enquire.....](#) [Read More.....](#)

**Description :**

Fits 0.2 mm diameter holes or larger, Straight



**Item Name :** BUTTON TIP MANIPULATOR

**Item Code :** 127-05

[Enquire.....](#) [Read More.....](#)

**Description :**

0.5 mm button tip, angled 10 mm  
May also be used for iris retraction



**Item Name :** IOL MANIPULATOR (FENZL)

**Item Code :** 127-06

[Enquire.....](#) [Read More.....](#)

**Description :**

Fits 0.2 mm diameter holes or larger, Angled 10 mm

