

SURGICAL KNIVES

Item Name : Sub 2.0 mm Phaco Slit Knives

Item Code : MM-1020

[Enquire.....](#) [Read More.....](#)

Description :

For Bi-manual phaco and micro-incision cataract procedures

MM-1014 : 1.4 mm, angled, full handle, bevel up

MM-1015 : 1.5 mm, angled, full handle, bevel up

MM-1016 : 1.6 mm, angled, full handle, bevel up

MM-1017 : 1.7 mm, angled, full handle, bevel up

MM-1018 : 1.8 mm, angled, full handle, bevel up

MM-1019 : 1.9 mm, angled, full handle, bevel up

MM-1020 : 2.0 mm, angled, full handle, bevel up



Item Name : Coaxial Microphaco Knife

Item Code : MM-1122

[Enquire.....](#) [Read More.....](#)

Description : MM-1122 : 2.2 mm, Angled, Full handle, double bevel



Item Name : Clear Comeal Knives, Parallel Sides

Item Code : MM-1332

[Enquire.....](#) [Read More.....](#)

Description :



MM-1327 : 2.75 Mm Angled, Full Handle, Double Bevel

MM-1330 : 3.0 Mm Angled, Full Handle, Double Bevel

MM-1332 : 3.2 Mm Angled, Full Handle, Double Bevel

Item Name : Stab Knives

Item Code : MM-2250

[Enquire.....](#) [Read More.....](#)

Description :

For side port or scleral incision

MM-1500 : 15⁰, Full handle

MM-2250 : 22.5⁰, Full handle

MM-3000 : 30⁰, Full handle

MM-4500 : 45⁰, Full handle



Item Name : Trapezoid Clear Comeal / Implant Knives

Item Code : MM-2832

[Enquire.....](#) [Read More.....](#)

Description :

MM-2832 : 2.85 mm - 3.2 mm

Angled, Full Handle, Bevel Up

MM-3032 : 3.0 mm - 3.2 mm Angled, Full Handle, Bevel Up

Dual-purpose design provides the ability to accomplish clear corneal and implant incisions with one keratome



Item Name : Blade Holders

Item Code : MM-4310



[Enquire.....](#) [Read More.....](#)

Description :

MM-4312 : 7.5cm Universal Handle

MM-4310 : 10cm Universal Handle

Item Name : Double Bevel Sideport Knives

Item Code : MM-5014

[Enquire.....](#) [Read More.....](#)



Description :

To create a self-sealing incision During bi-manual phaco procedures

MM-5010 : 1.0 mm, angled, full handle, double bevel

MM-5012 : 1.2 mm, angled, full handle, double bevel

MM-5014 : 1.4 mm, angled, full handle, double bevel

Item Name : Phaco Slit Knives

Item Code : MM-5510

[Enquire.....](#) [Read More.....](#)

Description :

For clear corneal or scleral incision

MM-5510 : 2.5mm, Angled, Full handle, Bevel up
MM-5568 : 2.65 mm, Angled, Full handle, Bevel up

MM-4738 : 2.75 mm, Angled, Full handle, Bevel up
MM-4740 : 2.85 mm, Angled, Full handle, Bevel up

MM-4741 : 2.80 mm, Straight, Full handle

MM-5570 : 3.0 mm, Angled, Full handle, Bevel up

MM-5571 : 3.0 mm, Straight, Full handle

MM-5572 : 3.2 mm, Angled, Full handle, Bevel up



MM-5573 : 3.2 mm, Straight, Full handle
 MM-5580 : 3.5 mm, Angled, Full handle,
 Bevel up

Item Name : MVR Knives

Item Code : MM-5563

[Enquire.....](#) [Read More.....](#)



Description :



For Vitreoretinal Infusion Port or
 Paracentesis Incision

MM-5555 : 19 Gauge, Full Handle

MM-5560 : 20 Gauge, Full Handle

MM-5563 : 23 Gauge, Full Handle

Item Name : MVR Knives Angled

Item Code : MM-5663

[Enquire.....](#) [Read More.....](#)



Description :



For Vitreoretinal Infusion Port or
 Paracentesis Incision

MM-5660 : 19 Gauge, Full Handle,
 Angled

MM-5661 : 20 Gauge, Full Handle,
 Angled

MM-5663 : 23 Gauge, Full Handle,
 Angled

Item Name : Needle Blades

Item Code : MM-5910

[Enquire.....](#) [Read More.....](#)



Description :

Versatile Blades useful in Pediatric
 Surgery, Capsulotomy, and Paracenteses

MM-5910 : 3.0mm

MM-5920 : 5.0mm

MM-5930 : 7.0mm

MM-5940 : 1.5mm

Item Name : Miniature Edged Blades

Item Code : MM-6200

[Enquire.....](#) [Read More.....](#)

Description :

MM-6200 : Chisel Tip

MM-6300 : Pointed Tip

MM-6400 : Round Tip

MM-6500 : Lance Tip

MM-6700 : Curved Tip

MM-6900 : Sharp all around

A variety of shapes and Edges for precise Microsurgical procedures.



Item Name : Lamellar Blades

Item Code : MM-6600

[Enquire.....](#) [Read More.....](#)

Description :

Flat, Polished Undersurfaces make ideal Undermining Blades

MM-6600 : 60°



Item Name : Spoon Knives

Item Code : MM-6821

[Enquire.....](#) [Read More.....](#)

Description :

For scleral tunnel incision

MM-6821 : 3.0 mm, angled, full handle, double bevel



Item Name : Full Handle Scleral Blades

Item Code : MM-6823H



[Enquire.....](#) [Read More.....](#)

Description :

MM-6823H

Item Name : Crescent Knives

Item Code : MM9402



[Enquire.....](#) [Read More.....](#)



Description : For scleral tunnel incision
MM-9402 : 12.8 mm, angled, Full handle,
bevel up Pack of 10 ETO Sterile