

Tonometer

Akriti

2011



Akriti Oculoplasty Logistics

**Email: sales@akriti.co.in,
oculoplasty.logistics@gmail.com**

Web : www.akriti.co.in

**1-284/OU/232 OU Colony, Shaikpet,
Hyderabad-500008, Andhra**

+91 40 3244 3781, 2356 9138 , 2356 9148

+91-40-2356 2729

TONOMETER

CODE: 30.010: - APPLANATION TONOMETER

APPLANATION TONOMETER TM-30.010 IS PRODUCED FOR THOSE WHO WISH THE TONOMETER TO REMAIN PERMANENTLY ON THE SLIT LAMP. IT IS MOUNTED ON A PIVOT ON THE MICROSCOPE AND FOR EXAMINATION IS SWUNG FORWARD IN FRONT OF THE MICROSCOPE.



CODE: 13109: - BERENS SCHIOTZ TONOMETER STERILIZER

OUR CONVENIENT SCHIOTZ TONOMETER STERILIZER CONSISTS OF A BASE WITH SLOTS TO STORE 7.5, 10 AND 15 GRAM WEIGHTS. THE REMOVABLE DISINFECTANT RESERVOIR CAP HAS A HOLE TO STORE THE TEST BLOCK AND ANOTHER HOLE THAT ALLOWS ONLY THE TONOMETER TIPS TO REST IN THE DISINFECTANT. THIS DESIGN PROTECTS THE REST OF THE TONOMETER FROM CONTACT WITH CLEANING AGENTS AND POSSIBLE CORROSION. THE STERILIZER IS EASY TO KEEP CLEAN AND DOUBLES AS A STORAGE CASE THAT ALLOWS FOR ALL TONOMETER PARTS TO BE KEPT TOGETHER AND READY TO USE. GULDEN'S STERILIZER IS MADE OF HIGH QUALITY, DURABLE PLASTIC THAT INSURES A LONG LIFE.



CODE: PT-H: - PERKINS HAND HELD APPLANATION TONOMETER

THE PERKINS HAS BEEN ESTABLISHED FOR WELL OVER 30 YEARS AS A 'WORLD STANDARD' FOR HAND HELD APPLANATION TONOMETERS.



CODE: TORC08-000: - TONOMETER REICHERT ATP NCT/PACH ST PROBE,

REICHERT'S NEW ATP AUTO NON-CONTACT TONOMETER / PACHYMETER IS THE FIRST FULLY INTEGRATED TONOMETER / PACHYMETER, COMBINING THE MEASUREMENT OF IOP AND CCT INTO ONE INSTRUMENT. STRIGHT TIP PROBE MODEL DESIGNED IN RESPONSE TO NUMEROUS CLINICAL STUDIES THAT HAVE SHOWN THE IMPORTANCE OF MEASURING CENTRAL CORNEAL THICKNESS AS WELL AS INTRAOCULAR PRESSURE IN THE MANAGEMENT OF GLAUCOMA PATIENTS, THE ATP INSTRUMENT COMBINES BOTH FUNCTIONS AND PROVIDES COMPLETE INFORMATION ON A SIMPLE-TO-READ, COLOR DISPLAY.



CODE: TOKE03-000: - KEELER PULSAIR EASYEYE TONOMETER,

KEELER'S PULSAIR EASYEYE IS THE ONLY HANDHELD NON CONTRACT TONOMETER WITH OVER 17 YEARS OF CLINICAL EXPERIENCE ON THE MARKET. THERE IS NO PRE-TEST PREPARATION, AND MINIMAL RISK OF CROSS INFECTION. THE PULSAIR EASYEYE ALSO PROVIDES EVEN BETTER CALIBRATION AND GREATER ACCURACY - ITS AUTOMATIC, PRECISION PERFORMANCE IS CALIBRATED TO READ WITHIN ± 1 MM HG OF GOLDMANN TONOMETRY.



PRODUCT DESCRIPTION

2011

CODE: R-900: - TONOMETER, HAAG-STREIT R-900 APPLANATION,

THE APPLANATION TONOMETER 900, AN ACCESSORY TO THE SLIT LAMP 900, CAN BE ESPECIALLY ADAPTED IN AN OPTIMUM MANNER TO COMBINE WITH OTHER HAAG-STREIT EQUIPMENT. BECAUSE OF ITS VERSATILITY, THE APPLANATION TONOMETER CAN ALSO BE MOUNTED ONTO AND COMBINED WITH OTHER EQUIPMENT.



CODE: TORC09-000: - TONOMETER, REICHERT TONOPEN AVIA

THE REICHERT TONO-PEN AVIA TONOMETER IS AN EASY-TO-USE, HAND-HELD INSTRUMENT THAT PROVIDES IOP READINGS THAT CORRELATE STRONGLY WITH GOLDMANN TONOMETRY. TONO-PEN AVIA HAS A LIGHTWEIGHT, ERGONOMIC DESIGN AND ADVANCED ELECTRONIC MEASUREMENT TECHNOLOGY THAT ENABLES OPERATORS TO TAKE FAST AND ACCURATE IOP MEASUREMENTS WITH MINIMAL TRAINING. TONO-PEN AVIA IS POWERED BY A LONG-LIFE TONO-PEN AVIA POWERCEL™ BATTERY PACK, GOOD FOR OVER 2,500 MEASUREMENTS. TONO-PEN AVIA™, LIKE THE TONO-PEN XL, USES SANITIZED OCU-FILM® TIP COVERS TO MINIMIZE THE RISK OF PATIENT CROSS CONTAMINATION.



CODE: TOME01-000: - REICHERT TONO-PEN® XL APPLANATION TONOMETER,

THE TONO-PEN XL APPLANATION TONOMETER IS AN EASY TO USE, HAND-HELD INSTRUMENT THAT PROVIDES IOP READINGS THAT CORRELATE CLOSELY WITH GOLDMANN TONOMETRY. TONO-PEN'S LIGHTWEIGHT ERGONOMIC DESIGN AND ADVANCED ELECTRONIC MEASUREMENT TECHNOLOGY ENABLE OPERATORS TO TAKE FAST AND ACCURATE IOP MEASUREMENT WITH MINIMAL TRAINING. WEIGHING ONLY 2.25 OUNCES, TONO-PEN XL IS EXTREMELY PORTABLE MAKING IT EASIER TO TAKE IOP MEASUREMENTS ANYWHERE ON ANY PATIENT.



CODE: 15169: - TIP CLEANER CADDY,

QUICKLY DISINFECTS TONOMETER TIPS REMOVABLE, REFILLABLE CUP EASILY MAINTAINED BY A TECHNICIAN SMALL FOOTPRINT OR WALL MOUNT KEEPS DROPS ORGANIZED NO NEED TO SEARCH FOR MISSING BOTTLES BUY ONE FOR EACH EXAM ROOM ONE SMALLER HOLE FOR HOLDING GONIOSOL UPSIDE DOWN PEROXIDE CUP HAS A FILL LINE TWO SMALLER HOLE FOR TONOMETER TIP DRYING

