

2011

Vision Charts



Akriti Oculoplasty Logistics


Hyderabad, India

1/1/2011

DISTANCE VISION CHARTS

CODE: 4636: - DIRECTIONAL 'E' CHART

THIS NEWLY REDESIGNED CHART IS MADE OF 11'X18' HIGH QUALITY 0.30' STYRENE. LINES INCLUDE 20/125, 20/100, 20/80, 20/60, 20/50, 20/40, 20/30, 20/25, 20/20, 20/15. MARKINGS PROVIDED FOR 20 FT EQUIVALENT SNELLEN, 6 METER AND LOGMAR. ROUNDED CORNERS.



CODE: 500010: - ETRS 10FT CHART

ETDRS CHART # 1 FOR 2000 CABINET, USE AT 10 FT (3M), 24.5X25.5 INCHES (62X65 CM), 20/200 TO 20/10.



CODE: 500011: - ETRS 10FT CHART # 2 FOR 2000 CABINET

ETDRS CHART # 2 FOR 2000 CABINET, USE AT 10 FT (3M), 24.5X25.5 INCHES (62X65 CM), 20/200 TO 20/10.



CODE: 500012: - ETRS 10FT CHART # 3 FOR 2000 CABINET

ETDRS CHART # 3 FOR 2000 CABINET, USE AT 10 FT (3M), 24.5X25.5 INCHES (62X65 CM), 20/200 TO 20/10.



CODE: 2121: - ETRS 4 METER CHARTS

TRANSLUCENT THESE CHARTS WORK GREAT EITHER IN THE CABINET OR HUNG ON A WALL! ETRS CHART # 1 RIGHT EYE, TRANSLUCENT



PRODUCT DESCRIPTION

2011

CODE: ESV3000: - ETDRS STANDARD ILLUMINATOR CABINET

ETDRS VIEWING CABINET, SELF-CALIBRATING, PHOTOPIC (85 CDM2) AND MESOPIC (3 CDM2) LIGHT TESTING, LED, AUTOMATIC 100 TO 240 VAC. WITH CUSTOMER PROVIDED PLUS (EXCEPT US). WITH WIRELESS REMOTE. CHART STORAGE IN BACK. MEASURES 24.5' H X 25.5' W X 3.75' D (62 W CM X 65 H CM X 36 CM D)



CODE: 4895: - FEINBLOOM DISTANCE LOW VISION FLIP CHART

FEINBLOOM DISTANCE LOW VISION FLIP CHART. USE AT 10 FT. COVERS 10/10 TO 10/700, ALL NUMBERS, LAMINATED PAPER WITH VINYL SPIRAL BOUND COVERS.



CODE: 800726: - HOTV LETTER CHART,

CHART IS SIMILAR TO THIS PHOTO




CODE: 600000: - ILLUMINATED CABINETS AND CHARTS

THE ILLUMINATED CABINET IS LIGHTWEIGHT AND VERSATILE AND USES 9'X14' TRANSLUCENT CHARTS. IT CAN BE PLACED ON A TABLE, HUNG ON A WALL OR PUT ON A STAND TO ACCOMMODATE ACCURATE 10 TO 20 FOOT VISUAL ACUITY SCREENING AND TESTING. HOUSED IN A DURABLE METAL CABINET, THE MODEL A PROVIDES DEPENDABLE SERVICE AND HAS BEEN THE MOST POPULAR CABINET OF PRACTITIONERS FOR OVER 25 YEARS. ONLY 5 POUNDS! USES F8T5D 8 WATT FLUORESCENT BULB. P/N 60600 ILLUMINATED CABINET DIFFERS SLIGHTLY IN THAT IT HAS A STORAGE COMPARTMENT FOR EXTRA CHARTS IN THE BACK.




CODE: 2205: - LANDOLT C FOLDING CHART

NOTATIONS FOR TESTING AT 3 METERS (10 FT.) ACUITY RANGE 30M TO 1.5M (20/200, 6/60 TO 20/10, 6/3 EQUIV.) TWO SMALLER GROUPS 7.5M TO 1.5M FOR REPETITIVE TESTING (20/50, 6/15 TO 20/10, 6/3 EQUIV.) MEASURES 45.5 X 45.5 CM (18 X 18 IN.) FOLDS TO 45.5 X 23 CM (18 X 9 IN.) CAN BE WALL-MOUNTED



CODE: 255000: - LEA CHART SET

LEA SYMBOL SET. 9X14 (23X35.6 CM) 20/100 TO 20/16, WITH RESPONSE KEY, FLASH CARDS, AND INSTRUCTIONS. USE AT 10 FT (3M).



CODE: 271100: - LEA NUMBER FOLDING 10FT DISTANCE CHART,

LEA NUMBERS FOLDING DISTANCE CHART FOR USE AT 10' (3M). PROVIDES 15 LINE TEST FOR YEAR OLDS AND UP. 18' X 20' (45.7 CM X 50.8 CM). FOLDS TO 18' X 10' (45.7 CM X 25.4 CM). HAS HOLE FOR HANGING.




CODE: 274000: - LEA NUMBERS SINGLE NUMBER BOOK FOR 10FT

LEA NUMBERS SINGLE NUMBER BOOK FOR USE AT 10 FEET (3 M). IDEAL FOR SITUATIONS IN WHICH THE CHILD/ADULT CANNOT PERFORM WHEN TESTED WITH A LINE TEST. THIS BOOK CAN HELP DETERMINE THE BEST RESOLUTION OF THE VISUAL PATHWAY IN A CHILD OR ADULT WITH AMBLYOPIA OR IMPAIRED VISION. INCLUDES ISOLATING MASK AND INSTRUCTIONS. 13 PAGES, 5' X 5' (12.7 CM X 12.7 CM) .




CODE: 255100 : - LEA SYMBOL 12 LINE TRANSLUCENT CHART SET,

LEA SYMBOL TEST SET. 12 LINE TRANSLUCENT DISTANT CHART SET. COMPLETE SET INCLUDES RESPONSE KEY, FLASH CARDS AND INSTRUCTIONS



CODE: 255001: - LEA SYMBOL CHART,

LEA SYMBOL CHART ONLY. 9X14 (23X35.6 CM) 20/100 TO 20/16, USE AT 10 FT (3M).



CODE: ESV1018: - LED NON-CALIBRATED CABINET WITH REMOTE,

LED NON-CALIBRATED PROFESSIONAL 10X18 INCH CABINET. AUTOMATIC 110 VAC. MEASURES 18.125' H X 10.125' W X 3' D, USES 10X18 INCH (25.4 X 45.7 CM) CHARTS. CONTROLLED BY WIRELESS REMOTE.



CODE: 259500: - MASSVAT 10 LINE LEA SYMBOL 10FT DISTANCE CHART,

FOR USE AT 10 FT, 20/40 AND 20/32 ONLY SIZE 3.5X9 SPIRAL - 4 PAGES, WITH RESPONSE KEY, FLASH CARDS



CODE: 259700 : - MASSVAT HOTV AND LEA SYMBOLS DISTANCE CHART

FOR USE AT 10 FT, 20/40 AND 20/32 ONLY, SIZE 10X10, CABINET, SET WITH RESPONSE CARD AND KEY


CODE: 259100: - MASSVAT LEA SYMBOL 10FT DISTANCE CHART

MASSVAT LEA SYMBOL DISTANCE CHART, FOR USE AT 10 FT, 20/80 TO 20/16, SIZE 9X14, FOR USE WITH SMALL ILLUMINATED CABINET, COMPLETE SET WITH RESPONSE KEY, FLASH CARDS, AND INSTRUCTIONS



CODE: 259111: - MASSVAT LEA SYMBOLS DISTANCE CHART SET

MASSVAT LEA SYMBOLS DISTANCE CHART, FOR USE AT 10 FT AND 20 FT, 20/80 TO 20/16, SIZE 9X14, FOR USE WITH SMALL ILLUMINATED CABINET, COMPLETE SET WITH RESPONSE KEY, FLASH CARDS AND INSTRUCTIONS.



CODE: 914000 : - PRE-SCHOOL REMOTE CONTROLLED ILLUMINATED CABINET,

PRE-SCHOOL REMOTE CONTROLLED ILLUMINATED CABINET, LED INSTA-LINE, WIRELESS CONTROLLED, RECHARGEABLE POWER PACK INCLUDED, 14.5H X 9.25W X 5D INCHES. (DOES NOT INCLUDE CHARTS, OCCLUDER OR HYPEROPIA GLASSES) 914000 REMOTE CONTROLLED ILLUMINATED CABINET:- REMOTE CONTROLLED ILLUMINATED CABINET, LED INSTA-LINE, WIRELESS CONTROLLED, RECHARGEABLE POWER PACK INCLUDED, INCLUDES LEA SYMBOLS, LEA NUMBERS, E, LETTER, AND HOTV CHARTS PLUS OCCLUDER, +1.75, +2.25 AND +2.50 HYPEROPIA GLASSES. 14.5H X 9.25W X 5D INCHES. 915000 REMOTE CONTROLLED ILLUMINATED CABINET, INSTA-LINE:- PRE-SCHOOL REMOTE CONTROLLED ILLUMINATED CABINET, LED INSTA-LINE, WIRELESS CONTROLLED, RECHARGEABLE POWER PACK INCLUDED, INCLUDES LEA SYMBOLS, AND HOTV CHARTS PLUS OCCLUDER, 14.5H X 9.25W X 5D INCHES. 916000 REMOTE CONTROLLED ILLUMINATED CABINET, LED INSTA-LINE:-



PRODUCT DESCRIPTION

2011

CODE: 800713: - SLOAN LETTER CHART DV,

DISTANCE CHART. SHOWS LETTERS D, V ON FIRST LINE. 10X18, TRANSLUCENT, WAS L10P




CODE: 800720: - SLOAN RANDOM LETTER CHART

DISTANCE CHART - RANDOM LETTERS 10X18 INCHES, TRANSLUCENT, WAS SL10P




CODE: 800714: - TUMBLING E DISTANCE CHART,

TUMBLING 'E' DISTANCE CHART, 10 FT, 10X18 INCHES. STYRENE, 20/16 TO 20/160. LINEAR SPACING



CODE: 600738: - TUMBLING 'E' CHART

TUMBLING 'E' CHART 10 FT (3M), PROPORTIONALLY SPACED, USABLE IN THE 9X14 INCH (23X36 CM) ILLUMINATED BOX OR HUNG ON A WALL.



MULTI-DISTANCE FOR PRESBYOPIC CORRECTIONS

PRODUCT DESCRIPTION

2011

CODE: PV4070-1: - COLENBRANDER READING, INTERMEDIATE AND COMPUTER CH,

CALIBRATED FOR 40 CM (16 IN), 63 CM (25 IN), 100 CM (40 IN)
MEASURES 8.5' X 11' INCLUDES COMPLETE INSTRUCTIONS AND RECORDING SHEET




CODE: PV-4000: - ENGLISH COLENBRANDER LOW VISION LETTER CHART,

FOR TESTING AT 1 METER (40 IN.) CHART CONTAINS SLOAN LETTERS IN LOGMAR INCREMENTS. OPTOTYPES RANGE IN SIZE FROM 50M TO 5M ON THE MAIN CHART AND FROM 6M TO 1M ON THE INSIDE CHART



CODE: PV4011: - ENGLISH COLENBRANDER LOW VISION READING CARD,


CARD CONTAINS PROPORTIONALLY SPACED READING SEGMENTS IN LOGMAR INCREMENTS RANGING IN SIZE FROM 10M TO 0.6M IN TIMES ROMAN FONT. FROM 4M ON, TWO SETS ARE PROVIDED TO AVOID MEMORIZATION. EACH LINE CONTAINS THE SAME NUMBER OF LETTERS AND SPACES



NEAR POINT CARDS


CODE: PV2106-90: - 90* SLOAN ETDRS FORMAT NEAR VISION CHART,

A 90 DEGREE TURN ON OUR MOST POPULAR NEAR VISION CHART! THIS FORMAT FACILITATES NATURAL LEFT TO RIGHT READING, WHILE LETTER SIZES DECREASE MAKING TESTING FASTER AND EASIER IN MANY CASES EACH SIDE IS DIFFERENT TO MINIMIZE MEMORIZATION



CODE: PV4041: - COLENBRANDER CONTINUOUS TEXT NEAR VISION CARD,

NOTATIONS FOR COMMONLY RECOMMENDED TESTING STANDARD OF 40 CM (16 IN.). ENGLISH READING CARD WITH SENTENCES OF EQUAL LENGTH AND DIFFICULTY LETTER SIZES FROM 6.3M TO .32M IN TIMES ROMAN FONT



PRODUCT DESCRIPTION


2011

NUMBERS PER LINE FOR 80, 60, AND 50 * ONE CHART
WITH 5 NUMBERS PER LINE FOR 40 AND 30 PLUS 6
NUMBERS FOR 25 AND 7 NUMBERS FOR 20, 15 AND 10

SINGLE CONTRAST LEVEL, VARIABLE SIZE>NEAR VISION 40CM (16 IN)

CODE: PV2118: - 10% SLOAN DOUBLE SIDED LOW CONTRAST VISION CHART

NOTATIONS FOR TESTING AT 40 CM (16 IN.) ACUITY
RANGE 8M TO .20M (20/400, 6/120 TO 20/10, 6/3 EQUIV.)
FRONT AND BACK SIDES ARE 10% LOW CONTRAST




CODE: PV4051: - COLENBRANDER CONTINUOUS TEXT CONTRAST CHART,

NOTATIONS FOR TESTING AT 40 CM (16 IN.) ACUITY
RANGE 8M TO .20M (20/400, 6/120 TO 20/10, 6/3 EQUIV.)
FRONT AND BACK SIDES ARE 10% LOW CONTRAST




CODE: PV2113-2B: - NEAR VISION POCKET CARD,

ACUITY RANGE 8M TO .20M (20/400, 6/120 TO 20/10, 6/3 EQUIV.).CALIBRATED FOR TESTING AT 40 CM (16 IN.)



CODE: PV2117: - SLOAN DOUBLE SIDED LOW CONTRAST VISION CARD,

NOTATIONS FOR TESTING AT 40 CM (16 IN.) ACUITY
RANGE 8M TO .20M (20/400, 6/120 TO 20/10, 6/3 EQUIV.)



CONTRAST ACUITY CHARTS:- SINGLE CONTRAST LEVEL, VARIABLE SIZE:- ANY-DISTANCE CHART

CODE: PC2129: - 1.25% CONTRAST SENSITIVITY CHART USING SLOAN OPTOT

FOR TESTING AT 3 METERS (10 FT.) ACUITY RANGE 25M TO 5M (20/160, 6/48 TO 20/10, 6/3 EQUIV.) TWO SMALLER GROUPS 6M TO 1.5M FOR REPETITIVE TESTING (20/40, 6/12 TO 20/10, 6/3 EQUIV.)



CODE: PV5000: - BAILEY LOVIE CHART SET

THE CHARTS STANDARDIZED THE TEST TASK SO THAT SIZE IS THE ONLY SIGNIFICANT VARIABLE FROM ONE ROW TO THE NEXT. THERE ARE 5 LETTERS PER ROW, BETWEEN-LETTER AND BETWEEN-ROW SPACINGS EQUAL TO LETTER SIZE, SIZE PROGRESSES BY A CONSTANT RATIO (0.1 LOG UNIT = 1,26X), AND LETTER DIFFICULTIES ARE BALANCED FOR EACH ROW.



CODE: PV4100: - COLENBRANDER MIXED CONTRAST CHART,

THE COLENBRANDER MIXED CONTRAST CHART CONTAINS BLACK AND 20% WEBER (11% MICHELSON) CONTRAST.




CODE: PV4061: - COLENBRANDER MIXED CONTRAST READING CARD,

DESIGNED AS A SIMPLE SCREEN TEST. THE USUAL QUESTION FOR A READING TEST IS: 'WHAT IS THE SMALLEST (BLACK) LINE YOU CAN READ?' ON THE CONTRAST CARD ONE ADDS THE QUESTION 'WHAT IS THE SMALLEST (GREY) LINE YOU CAN READ?' THE DIFFERENCE BETWEEN THE TWO MEASUREMENTS WILL BE OBVIOUS AND CAN SIMPLY BE EXPRESSED IN 'LINES'. PRELIMINARY FINDINGS INDICATE A 2 OR 3 LINE DIFFERENCE IS NORMAL.



CODE: PV2057: - HOTV LOW CONTRAST FLIP BOOK,

MULTI-DISTANCE TESTING CONTAINS FIVE PAGES OF OPTOTYPES OF 9.5M SIZE WITH 5 SLOAN LETTERS CONTRAST VARIES BY PAGE: 25%, 10%, 5%, 2.5%, AND 1.2%




CODE: PV2180: - MULTI-CONTRAST SLOAN LETTER FLIP BOOK,

FOR TESTING AT 2 METERS (6.5 FT.) EACH PAGE IT'S OWN CONTRAST LEVEL AT 100%, 2.5% AND 1.25% ACUITY RANGE 20M TO 1.6M




CODE: PV2511: - PATTI PICS™ LOW CONTRAST FLIP BOOK,

MASSVAT DESIGNED BY: LUISA MAYER, PH.D. AND BRUCE MOORE, O.D. MULTI-DISTANCE TESTING EACH FLIP CHART PAGE HAS FIVE SYMBOLS IN THE MASSVAT FORMAT AT 9.5M SIZE THE CONTRAST LEVELS VARY BY PAGE: 25%, 10%, 5%, 2.5% AND 1.2%




CODE: PV2512: - PATTI PICS™ LOW CONTRAST FLIP BOOK SET 2,

MASSVAT DESIGNED BY: LUISA MAYER, PH.D. AND BRUCE MOORE, O.D. MULTI-DISTANCE TESTING EACH FLIP CHART PAGE HAS FIVE SYMBOLS IN THE MASSVAT FORMAT AT 9.5M SIZE THE CONTRAST LEVELS VARY BY PAGE: 25%, 10%, 5%, 2.5% AND 1.2%




CODE: PV2705: - PV NUMBERS™ LOW CONTRAST FLIP BOOK

MULTI-DISTANCE TESTING EACH FLIP CHART HAS PAGES OF FIVE NUMBERS IN MASSVAT FORMAT AT 9.5M SIZE THE CONTRAST LEVELS VARY BY PAGE: 25%, 10%, 5%, 2.5%, AND 1.2%.



CODE: PV2704: - PV NUMBERS™ LOW CONTRAST FLIP BOOK 1

MULTI-DISTANCE TESTING EACH FLIP CHART HAS PAGES OF FIVE NUMBERS IN MASSVAT FORMAT AT 9.5M SIZE THE CONTRAST LEVELS VARY BY PAGE: 25%, 10%, 5%, 2.5%, AND 1.2%.



PRODUCT DESCRIPTION

2011

**CODE: PV2192: - SLOAN LOW CONTRAST FLIP BOOK,
PRICE: 0.00**

SCRAMBLED VERSION OF CAT.NO. 2190 MULTI-DISTANCE TESTING CONTAINS FIVE PAGES OF OPTOTYPES OF 9.5M SIZE WITH 5 SLOAN LETTERS CONTRAST VARIES BY PAGE: 25%, 10%, 5%, 2.5%, AND 1.2%

CODE: PV2191: - SLOAN LOW CONTRAST FLIP BOOK

MULTI-DISTANCE TESTING CONTAINS FIVE PAGES OF OPTOTYPES OF 9.5M SIZE 5 SLOAN LETTERS CONTRAST VARIES BY PAGE: 25%, 10%, 5%, 2.5%, AND 1.2%

**CODE: PV2125: - SLOAN LOW CONTRAST LETTER SET
BOOK**

NOTATIONS FOR TESTING AT 1 METER (40 IN.) ACUITY RANGE 20M TO 1.6M (20/200, 6/120 TO 20/32, 6/9.5 EQUIV.) CONTRAST VARIES BY PAGE: FULL CONTRAST, 25% , 10%, 5%, 2.5%, 1.25%, AND .06% SPIRAL BOUND



**CODE: PV2380: - TUMBLING 'E' LOW CONTRAST FLIP
BOOK,**

FOR TESTING AT 2 METERS (6.5 FT.) EACH PAGE IT'S OWN CONTRAST LEVEL AT 100%, 2.5% AND 1.25% ACUITY RANGE 20M TO 1.6M

CONTRAST ACUITY CHARTS:- SINGLE CONTRAST LEVEL, VARIABLE SIZE > 4 METER

**CODE: PV2212: - LANDOLT 'C' TRANSLUCENT LOW CONTRAST
CHART 10%,**

NOTATIONS FOR TESTING AT 4 METERS (13 FT.). ACUITY RANGE 40M TO 2M (20/200, 6/60 TO 20/10, 6/3 EQUIV.). MEASURES 62 X 65 CM (24.5 X 25.5 IN.)



PRODUCT DESCRIPTION

2011

CODE: PV2023: - MASSVAT HOTV LOGARITHMIC VISUAL ACUITY CHART, PRICE: 0.00, SHIPING COST: 0.00

FOR TESTING AT 40 CM (16 IN.) MEASURES 9.5 X 16 CM (3.75 X 6.25 IN.) PROTECTIVE COATING ALLOWS WASHING OR CLEANING WITH ALCOHOL WITHOUT ANY DAMAGE TO CHART FACES.




CODE: PV2113-2: - NEAR VISION POCKET CARD, PRICE: 0.00, SHIPING COST: 0.00

ACUITY RANGE 8M TO .20M (20/400, 6/120 TO 20/10, 6/3 EQUIV.).




CODE: PV2594: - PATTI PICS RUNGE POCKET NEAR VISION TEST CARD, PRICE: 0.00, SHIPING COST: 0.00

ACUITY RANGE 8M TO .20M (20/400, 6/120 TO 20/10, 6/3 EQUIV.) MEASURES 9.5 X 16 CM (3.75 X 6.25 IN.)




CODE: PV2586: - PATTI PICS™ MASSVAT® NEAR VISION POCKET CARD, PRICE: 0.00, SHIPING COST: 0.00

MASSVAT® DESIGNED BY: LUISA MAYER, PH.D. AND BRUCE MOORE, O.D. NOTATIONS FOR VISUAL ACUITY TESTING AT 40 CM (16 IN.) MASSVAT FORMAT BOTH SIDES OF NEAR POINT VISION TESTING CARD. CLEAR COATING TO PROTECT EYE CHARTS, WILL NOT BECOME DAMAGED IF WASHED OR CLEANED WITH ALCOHOL.



CODE: PV2509: - PATTI PICS™ POCKET SIZE NEAR VISION CARD, PRICE: 0.00, SHIPING COST: 0.00

MASSVAT DESIGNED BY: LUISA MAYER, PH.D. AND BRUCE MOORE, O.D. NOTATIONS FOR TESTING AT 40 CM (16 IN.) MEASURES 9.5 X 16 CM (3.75 X 6.25 IN.)



CODE: PV2748: - PV NUMBERS RUNGE POCKET NEAR VISION TEST CARD, PRICE: 0.00, SHIPING COST: 0.00

ACUITY RANGE 8M TO .20M (20/400, 6/120 TO 20/10, 6/3 EQUIV.) MEASURES 9.5 X 16 CM (3.75 X 6.25 IN.) DEVELOPED BY DR. PETER RUNGE, MD.

